

Chairside machines

# What are the differences?



CORiTEC one



CORiTEC one+

Number of axes	4	5
Number of tool positions	6	10
Spindle power / rpm	750 watt / 100.000 rpm	
Quick clamping system (maually)	✓	✓
Dry processing	✓*	✓*
Wet processing	✓	✓
98 mm blank processing	-	✓
CAD/CAM blocks / glass ceramics	✓	✓
Premilled abutments	✓	✓
Control-PC	integrated	integrated
10'1" HD touchscreen	✓	✓
Integrated coolant tank	integrated	integrated
LED statusbar	✓	✓
Auto calibration set	✓**	✓**
Compressed air	no external air required	no external air required
Collet	3 mm stainless steel	3 mm stainless steel
C-Clamp processing (optionally)	-	✓
Integration in exocad® chairside workflow	✓	✓
Turning processing (universal fixture required)	-	✓

Chairside machines  
**Adapters**



**CORiTEC one**



**CORiTEC one+**

98 mm blank holder	-	included
Universal fixture for glass ceramic and premilled abutments	included	included
6-fold glass ceramic holder (for 98 mm blank holder)	-	included
C-Clamp blank holder	-	available as accessory
Nt-trading 6-fold holder premilled abutment inkl. change bars	-	available as accessory
Medentika 6-fold holder premilled abutment inkl. change bars	-	available as accessory
DESS 6-fold holder premilled abutment inkl. change bars	-	available as accessory
Zfx 6-fold holder premilled abutment inkl. change bars	-	available as accessory
Glass ceramic adapter for universal fixture	included	included
Nt-trading fixture for universal adapter	available as accessory	available as accessory
Medentika fixture for universal adapter	available as accessory	available as accessory
DESS fixture for universal adapter	available as accessory	available as accessory
Zfx fixture for universal adapter	available as accessory	available as accessory